


CLAVE[®]

INTRAVESICAL ADMINISTRATION DEVICE

For the Safe Delivery of Hazardous Medications


Now you can have the ultimate in flexibility when conducting Intravesical Chemotherapy Administration. The CLAVE Connector offers a fail safe, closed connection device that will prevent spills, leaks and exposure.

CLAVE[®] Connector is Mechanically and Microbiologically Closed



Graduated Adapter attaches to any urinary catheter

 **Safe Handling Made Simple**



Use the Spiros[™] Closed Male Connector on the infusion device to complete the closed system



ICU Medical, Inc.

951 Calle Amanecer, San Clemente, CA 92673, U.S.A.
Tel: +1 949 366 2183 • www.icumed.com

US and Foreign Patents Pending. M1-1216 Rev. 1

Intravesical Administration Device



 **Safe Handling Made Simple**