



SwabTip™

Male-Luer Disinfecting Cap

Enhance patient safety with the only male-luer disinfecting cap that covers the entire open male luer.

icumedical
human connections

Provide Closure to Help in Infection Control



SwabTip™

Male-Luer Disinfecting Cap

SwabTip is a male-luer disinfecting cap that covers the entire male luer, disinfecting the luer tip with 70% IPA while providing a physical barrier to protect the internal luer threads. SwabTip works fast, disinfecting the open male luer in just 30 seconds and providing continuous protection for up to 4 days if not removed.¹

Dispensed from a convenient foil strip that can easily hang from any IV pole, you can get easy access to male luer disinfection where and when you need it.

To learn how SwabTip can help protect your patients,
call 866.829.9025 or visit www.icumed.com

Ordering Information

Model #	Description	Strips/Case	Strips/Box	Caps/Strip
SWTP-5	Disinfecting Cap for Male Luers	200	50	5

icumedical
human connections

¹ ICU Medical 510(k) Summary K190918