



# Deltaven<sup>®</sup> Fast Flash<sup>™</sup>

Closed IV Catheter System

Integrated peripheral IV catheter system offering closed system blood control and passive safety needle safety technology.

**icumedical**  
human connections

# Help Minimize Blood Exposure While Reducing the Risk of Accidental Needlestick Injury

The Deltaven Fast Flash Closed IV Catheter System combines closed system blood control and passive safety needle technology to help your efforts to minimize occupational exposures and comply with Occupational Safety and Health Administration (OSHA) bloodborne pathogen standards.<sup>1</sup>



## Help minimize the risk of blood exposure

Avoid blood leakage during catheter insertion and length of dwell with integrated closed system blood control.



## Reduce the risk of needlestick injury

Avoid accidental needlesticks and needle exposure with passive safety needle technology.



## Avoid multiple insertion attempts

Help confirm vessel entry with notched needle technology and dual vented system.

## Product Information

Single Port with End Caps	Dual Port with End Caps	Dual Port with Needlefree Connector	Gauge	Length	Gravity Flow (mL/min)	HPI Flow (mL/sec) <sup>2</sup>	HPI Rating (psi)
3804773	3800573	3800473	● 26G	0.75 in	15	–	–
3814773	3811073	3811473	● 24G	0.75 in	24	4	325
3825773	3821573	3821473	● 22G	1.00 in	31	5	325
3824773	3820573	3820473	● 22G	0.75 in	37	5	325
3837773	3833573	3833473	● 20G	1.75 in	51	7	325
3836773	3832573	3832473	● 20G	1.25 in	56	8	325
3835773	3831573	3831473	● 20G	1.00 in	59	8	325
3847773	3843573	3843473	● 18G	1.75 in	78	>10	325
3846773	3842573	3842473	● 18G	1.25 in	83	>10	325
3845773	3841573	3841473	● 18G	1.00 in	86	>10	325
3856773	3852573	3852473	● 16G	1.25 in	125	>10	325

Deltaven Fast Flash is a registered trademark of Delta Med SpA

<sup>1</sup> United States Occupational Safety and Health Administration (OSHA). Bloodborne Pathogens Standards. 29 CFR 1910.1030.

**icumedical**  
human connections