



# Gripper Micro™

## Port Access Needle

### Blunt Cannula, Safety Needle

- › A safety port access needle to administer or withdraw fluids from implanted ports
- › Available in different sizes and lengths
- › Non-sharp, blunt cannula for safe removal and disposal
- › Trocar, non-coring needle that locks into place for sharps disposal with an audible and visible click
- › Tiny footprint with thick cushioned needle platform sits against the skin when in-situ
- › Non-siliconized needle
- › The 20 cm-long extension set is available with or without Y-site with luer-activated needleless injection site
- › Extension set with color-coded clamp(s)
- › Priming volume: 0.3 mL or 0.8 mL depending on the extension set with or without the Y-site versions
- › The blunt cannula design contributes to faster flow rate
- › Pressure-rated up to 45 psi
- › Positive ANTT evaluation

### Ordering Information

#### Needleless Y-site — luer-activated needleless injection site

Product Code	Gauge	Length	Units
21-3271-24		¾ in	
21-3272-24	● 22	1 in	12/box
21-3273-24		1 ¼ in	
21-3276-24		¾ in	
21-3277-24	● 20	1 in	12/box
21-3278-24		1 ¼ in	
21-3280-24		¾ in	
21-3281-24	● 19	1 in	12/box
21-3282-24		1 ¼ in	

#### Without Y-site

Product Code	Gauge	Length	Units
21-3251-24		¾ in	
21-3252-24	● 22	1 in	12/box
21-3253-24		1 ¼ in	
21-3256-24		¾ in	
21-3257-24	● 20	1 in	12/box
21-3258-24		1 ¼ in	
21-3260-24		¾ in	
21-3261-24	● 19	1 in	12/box
21-3262-24		1 ¼ in	

